

8400B AUTOSAL LABORATORY SALINOMETER

The Best in Precision Salinometer Measurements



GUILDLINE INSTRUMENTS 8400B HAS BEEN THE LEADING SALINITY MEASUREMENT INSTRUMENT FOR OVER 20 YEARS. This Salinometer employs a unique continuous flow system, where the sample water is drawn under low air pressure from the original sample bottle. A high stability temperature control bath and heat exchanger maintain the sample at a precisely defined temperature during analysis, avoiding the need for temperature compensation. The accuracy of better than 0.002 Equivalent Practical Salinity Units (PSU) has been achieved by improvements in the temperature stability of the bath.

The unique continuous sample flow system offers a marked improvement in speed and convenience over the conventional salinometer method requiring several discrete sample volumes for flushing and measurement.

Only 50 milliliters of sample water are required to achieve a final reading, including flushing. Semi-automatic operation gives the final digital readout of conductivity ratio in less than one minute, with a resolution of better than 0.0002 Equivalent PSU.

The 8400B has a redesigned sample bottle holder, which provides for easy attachment of original sample bottles. There is a bottle data logging facility built in as standard, which provides output of conductivity ratio, bottle number, bath temperature and standardization reference number.

FEATURES

- Accuracy < 0.002 Salinity
- Range 2 to 42 in Salinity
- Improved Temperature Stability
- Small Samples and Rapid Measurements
- Data Logging Standard
- Improved Maintainability

The 8400B "Autosal" provides incredibly accurate salinity measurements in under a minute using just 50 ml of sample water

Bath temperature is selectable from 18 °C to 33 °C in 3 °C steps. There is a window for viewing sample water under test flowing through the conductivity cell. Operating controls include conductivity range select, display function switch for zero, standby or read, bath temperature control check, standard seawater standardization set, temperature bath drain and fill controls, and a sample flow rate control which permits deep cold water samples to be measured immediately.

The 8400B has many improvements in the design and layout of the instrument. Particular care has been taken to improve the power supply with better input filtering and a new cable harness. There is a heater failure lamp included on the front panel, and a front sub-panel reduces wear on the internal hoses, when the front door is opened for maintenance. There are also significant changes internally to improve ease of maintenance, such as access to the pumps.

8400B General Specifications

| | |
|------------------------------|---|
| Measurement Range: | 0.0001:1.15 Conductivity Ratio 0.004 to 76mS/cm 2 to 42 Equivalent Practical Salinity Units (PSU) |
| Accuracy: | < ± 0.0001 Conductivity Ratio, @ same set point temperature as standardization and within -2 °C and +4 °C of ambient. By calculation & substitution in the Bennett equation or the UNESCO tables, < ± 0.002 Equivalent PSU. |
| Short Term Stability: | < ± 0.00005 for 24 hours without re-standardization |
| Maximum Resolution: | < 0.00001 Conductivity Ratio < 0.0002 mS/cm @ 15 °C and 35 PSU < 0.0002 Equivalent PSU |
| Sample Volume: | Maximum required – 100 milliliters (starting from fresh water in the cell), including flushing volume. About 50 milliliters for 3 Equivalent PSU difference in samples. |
| Scale Suppression: | Linear scale of conductivity ratio having 22 steps from 0 to 2.2, where 2.0 corresponds to seawater of 35 Equivalent PSU. Maximum reading is 2.29999, corresponding to approximately 42 Equivalent PSU. |
| Bath Temperature: | Selectable from 18 °C to 33 °C in 3 °C steps, accuracy ±0.02 °C, stability ± 0.001°C per day. Selected temperature should be within (ambient +4) °C and (ambient -2) °C. |
| Outputs: | TTL compatible BCD outputs of numerical display readings and thumbwheel switch bottle logger. |
| Maintenance: | Cell easily cleanable by bottle brush. Can be removed for cleaning and replaced without change in calibration. |
| Water Bath Volume: | 4.4 U.S. gallons (16.8 liters). |
| Power Required: | 115 or 230 volts ± 10%, 50/60 Hz, 400 watts maximum |
| Temperature: | Operating: 16 °C to 37 °C Storage: -40 °C to 70 °C |
| Dimensions: | H 678 mm (26.7 in); W 533 mm (21 in); D 559 mm (22 in) |
| Weight: | Bath Empty: 50 kg (110 lbs) Bath Full: 70 kg (150 lbs) |

ORDERING INFORMATION

| | |
|---------------------|----------------------------------|
| 8400B | "Autosal" Laboratory Salinometer |
| Accessories: | |
| 84001 | Conductivity Cell |
| 84004 | Thermistor Kit |
| 84006 | Pump Motor Replacement Kit |
| 84007 | Maintenance Kit |
| 84105 | External Sample Pump |

GUILDLINE IS DISTRIBUTED BY:

Guildline Instruments Limited
Smiths Falls, Ontario, Canada, K7A 4S9
Phone: (613) 283-3000 Fax: (613) 283-6082
Web: www.guildline.com

18970-01-85 Rev. B1 Copyright © 2022.02.01 Guildline Instruments Limited. All rights reserved. Subject to change without notice.